






INSULIN PENS



<p>AUTOPEN CLASSIC <i>Oceen Mumford</i> Reusable</p> 	<p>HUMALOG KWIKPEN <i>Eli Lilly</i> Disposable</p> 	<p>HUMAPEN LUXURA HD <i>Eli Lilly</i> Reusable</p> 	<p>HUMAPEN MEMOIR <i>Eli Lilly</i> Reusable</p> 	<p>ORIGINAL PREFILLED PEN <i>Eli Lilly</i> Disposable</p> 	<p>FLEXPEN <i>Novo Nordisk</i> Disposable</p> 	<p>NOVOPEN 3 <i>Novo Nordisk</i> Reusable</p> 	<p>NOVOPEN JUNIOR <i>Novo Nordisk</i> Reusable</p> 	<p>SOLOSTAR <i>Sanofi-Aventis</i> Disposable</p> 
<p>INSULIN USED Humalog</p> <p>INSULIN DELIVERY Comes in a 1-unit increment version, which can dispense up to 21 units at a time, and a 2-unit increment version, which can dispense up to 42 units at a time</p>	<p>INSULIN USED Humalog, Humalog Mix 75/25, Humalog Mix 50/50</p> <p>INSULIN DELIVERY One-unit increments, up to 60 units at a time</p> <p>DETAILS Pen is blue, with labels of different colors for various types of insulin.</p>	<p>INSULIN USED Humalog</p> <p>INSULIN DELIVERY Half-unit increments, up to 30 units at a time</p> <p>DETAILS Can dispense half units of insulin (suitable for children). Pen is green and made out of metal.</p>	<p>INSULIN USED Humalog</p> <p>INSULIN DELIVERY One-unit increments, up to 60 units at a time</p> <p>DETAILS Has a digital display with memory, which allows users to see the time, date, and amount of the last 16 doses delivered. Nonreplaceable battery has a three-year lifespan. Pen is burgundy and made out of metal.</p>	<p>INSULIN USED Humulin N, Humulin 70/30</p> <p>INSULIN DELIVERY One-unit increments, up to 60 units at a time</p> <p>DETAILS Pen is white, with labels of different colors for various types of insulin. (This pen is no longer sold for use with Humalog, Humalog Mix 75/25, and Humalog Mix 50/50.)</p>	<p>INSULIN USED Levemir, NovoLog, NovoLog Mix 70/30</p> <p>INSULIN DELIVERY One-unit increments, up to 60 units at a time</p> <p>DETAILS Pen is blue, with labels of different colors for various types of insulin.</p>	<p>INSULIN USED NovoLog</p> <p>INSULIN DELIVERY One-unit increments, from 2 to 70 units at a time</p> <p>DETAILS Available in silver, green, or blue. Pen is made out of metal.</p>	<p>INSULIN USED NovoLog</p> <p>INSULIN DELIVERY Half-unit increments, up to 35 units at a time</p> <p>DETAILS Can dispense half units of insulin (suitable for children). Pen is blue with colorful detailing.</p>	<p>INSULIN USED Apidra or Lantus</p> <p>INSULIN DELIVERY One-unit increments, up to 80 units at a time</p> <p>DETAILS Apidra and Lantus versions of this pen have different colors (blue for Apidra, gray for Lantus) and feel different to the touch to help users distinguish between the types of insulin.</p>

DELIVERY DETAILS
Delivery button is on the side rather than on the end. Insulin is delivered with a spring mechanism, so the pressure required to push the button is the same regardless of dose size or needle gauge. Pen is available in blue and green.